



Novanta
PHOTONICS

WHITEPAPER

LASER TWEEZERS RAMAN SPECTROSCOPY (LTRS)

Reduce the Operator Input Needed to Acquire Single Cell Raman Spectra

Laser Tweezers Raman Spectroscopy (LTRS)

Reduce the Operator Input Needed to Acquire Single Cell Raman Spectra

Abstract

One in two people born after 1960 in the UK will be diagnosed with some form of cancer during their lifetime. Analysis of cancer tissue at the bulk level means individual cell attributes are averaged, and therefore single cell methods are of interest, both in a diagnostic setting and in the context of drug-cell interaction studies. In recent years, significant progress has been made into the label-free detection and discrimination of individual cancer cells using Laser Tweezers Raman Spectroscopy (LTRS). However, LTRS methods have invariably involved a high degree of manual intervention, and before this technique can be translated into a clinical setting a greater degree of automation is required.

Introduction

First observed in 1928 by Raman and Krishnan, Raman scattering (RS) is the inelastic scattering of light from a sample, which reveals information relating to the distribution of vibrational states (Figure 1). The non-destructive, label-free nature of RS makes the technique particularly attractive for the analysis of biological samples. However, the low efficiency at which Raman scattering occurs compared to fluorescence is a fundamental limitation, making the choice of substrate and sample preparation paramount. With typical signal collection times in the order of seconds, a method of immobilizing samples during the acquisition of a spectrum is required, and the majority of biological Raman studies involve single cells either cultured or spun onto substrates.




Figure 1: Energy level diagram showing vibrational states probed by Raman spectroscopy

An alternative method of maintaining the position of micron-sized objects including individual cells is optical trapping, first demonstrated by Ashkin in 1986. Optical trapping of cells has been achieved with both a single beam (laser tweezers, LT) and using two counter propagating beams, in both cases requiring the generation of a high intensity optical gradient, resulting in forces of the order of pico Newtons. For a single laser beam focused by a high numerical aperture lens, these forces can be sufficient to hold an individual cell above a substrate. As such, LT provide two key advantages for single cell analysis: first, cells can be isolated away from a Raman active/fluorescent substrate, and second, it is possible to

manipulate a trapped object, allowing cells of interest to be isolated from a larger population.

Laser Tweezers Raman Spectroscopy (LTRS)

The combination of LT and RS in a single experimental arrangement is straightforward, since the low signal levels found in RS generally require the use of a microscope objective and the means to couple laser radiation to a sample. The LTRS arrangement outlined utilizes a ventus 1064 together with a ventus 532 for optical trapping and Raman excitation respectively. The use of a low cost microfluidic flow chamber controlled via a syringe pump allows individual cells to be automatically trapped while a Raman spectrum is acquired, and greatly reduces the amount of operator input required compared to traditional methods of optical trapping.

Circulating Tumour Cells (CTCs) are cells which are shed into the blood system from a primary tumour. In the case of prostate cancer, prognosis is strongly dependent on whether cancer has spread to secondary sites, with outlook particularly poor for individuals with secondary bone metastases. It has been suggested that, for late stage (androgen independent) prostate cancer, detecting the number of CTCs in a sample of peripheral blood may also have some diagnostic potential, both in terms of establishing the current staging of the disease and in assessing the effectiveness of treatment. To verify the potential for LTRS to be used to differentiate between CTCs and white blood cells, the metastatic PC-3 prostate cell line and Jurkat cells were analyzed using the LTRS system outlined in Figure 2.




Figure 2: The experimental set-up

As shown in Figure 3, the average PC-3 Raman spectrum was found to display increased intensities representative of cytochrome-c in addition to enhanced protein/lipid bands. The Jurkat cells were found to display increased intensities indicative of nucleic acids. Using Linear Discriminate Analysis, classification rates in excess of 95% were obtained for the two cell lines. Furthermore, using only the ratio of Raman bands representative of cytochrome-c and nucleic acids, the two cells lines could be successfully classified at acquisition times as short as 2 seconds.




Figure 3: The average PC-3 Raman spectrum, showing increased intensities representative of cytochrome-c in addition to enhanced protein/lipid bands

Following the classification of PC-3 and Jurkat cells, a fully automated method of trapping and Raman analysis was explored using a second microfluidic flow chamber consisting of two individual channels. By alternating the position of the microfluidic flow cell relative to the focused beams, time-dependent changes in levels of cytochrome-c were identified for a single cell line exposed to different conditions.

LTRS arrangements can utilize a single laser for both trapping and Raman excitation, or two separate wavelengths. Trapping of individual cancer cells was achieved using Laser Quantum's ventus 1064 laser. The ventus 1064 outputs a Single Transverse Mode (TEM00) with an M2 value of <1.4 as required for the generation of an intense optical gradient. The optically trapped cells are then illuminated by Laser Quantum's ventus 532 laser. Illumination at 532 nm results in enhanced (resonant) Raman scattering by cytochrome-c, a heme protein involved in cellular respiration, allowing the two cell lines to be distinguished.

An important factor for Raman spectroscopy is the line width of the excitation laser, which will impact on how well closely-spaced Raman bands may be resolved. Typically, Raman line widths within single-cell spectra are of the order of 10 cm⁻¹ and as a result the 30 GHz line width of the ventus 532 (equal to 1 cm⁻¹) is ideal. For any dual wavelength LTRS systems, a further consideration is the pointing stability of the two beams. Both the ventus 1064 and ventus 532 have excellent pointing stability (<10 μrad/°C), ensuring the two beams remain focused at the same position over extended time frames.

The combination of a low cost microfluidic flow chamber and dual wavelength LTRS system has greatly reduced the level of operator input needed to acquire single cell Raman spectra, whilst supplying invaluable data for cancer research. The experiments outlined reveal the potential for moving toward a single automated system, requiring simultaneous control of spectrometer, syringe pump, microscope translation stage and lasers. Both the ventus 1064 and ventus 532 lasers can be controlled via an RS232 interface or using LabVIEW drivers, allowing lasers to be fully integrated into a single experimental set-up.

Novanta Benefits

Novanta is uniquely positioned to solve even the most complex challenges for OEMs, system integrators, and end-use customers seeking to advance their manufacturing processes with high precision laser systems. With some of the most well-known brands in the industry and in-country application and service support, Novanta delivers reliable, precise, and durable components and sub-systems.

Our Applications Testing Labs offer application and proof-of-concept testing to OEMs, system integrators, material manufacturers, processors, and end-users of automated machinery. Novanta Application Engineers are laser processing experts, and understand the parameters that will ensure successful, efficient laser processing. Using laser and beam steering equipment from well-known Novanta brands, our Application Engineers will determine the key product parameters and processing know-how to achieve the desired results.

Acknowledgments

Thank you to Dr. S. Casabella, Dr. P. Scully, Prof. N. Goddard, Prof. P. Gardner from The University of Manchester
Reference: S. Casabella et. al. Automated analysis of single cells using Laser Tweezers Raman Spectroscopy. Analyst, 2016, 141, 689

CONTACT US

E photonics@novanta.com | photonics.china@novanta.com | photonics.japan@novanta.com

W novantaphotonics.com